

CareNavigatorPH

Navigate to the right care, at the right time.

A mobile healthcare navigation app for Filipino communities

SDG 3 | Good Health and Well-Being

Student Name [Last Name], [First Name]

Course / Section Advanced Mobile Programming / [Section]

Instructor [Instructor Name]

Date _____

Filename Week1_Activity_LastName_FirstName

PROJECT SNAPSHOT

CareNavigatorPH helps people describe symptoms safely, understand the appropriate level of care, find verified hospitals with relevant services and availability, and begin a consultation without presenting AI guidance as a medical diagnosis.

Advanced Mobile Programming | CareNavigatorPH

1. Community Problem & SDG Alignment

PROBLEM STATEMENT

People in many Philippine communities struggle to decide where to seek care because hospital information is fragmented, service availability is unclear, and symptoms can be difficult to interpret. This can cause delayed treatment, unnecessary travel, crowded facilities, and confusion during urgent situations.

Why this matters in the community

- Patients may not know which nearby facility has the right department, specialist, emergency capability, beds, or rooms.
- Rural residents, older adults, caregivers, and first-time patients may spend extra time calling or visiting multiple facilities.
- Online health information can be misleading when it lacks local context, safety warnings, and a clear path to professional care.

Selected UN Sustainable Development Goal

GOOD HEALTH AND WELL-BEING

Ensure healthy lives and promote well-being for all at all ages.

SDG

- 3 **Direct alignment:** Target 3.8 focuses on access to quality essential healthcare services. CareNavigatorPH improves access by connecting users to appropriate, verified, and locally available care.

Intended impact

The app aims to shorten the path from uncertainty to appropriate professional care. It supports informed decisions and timely referrals while keeping emergency escalation, privacy, and clinician review central to the experience.

Problem and SDG Alignment | 2

2. Project Proposal

Core concept

CareNavigatorPH is a mobile healthcare navigation platform designed for Philippine communities. A user can describe symptoms, receive preliminary guidance about urgency and the type of care to seek, compare verified hospitals, and start a consultation. The app does not diagnose or replace a licensed healthcare professional.

Target audience

PRIMARY USERS

- Patients and caregivers
- Residents searching for nearby care
- People unsure which department or hospital level they need

SECONDARY USERS

- Doctors and hospital staff
- Hospital administrators
- Community health workers and referral partners

Primary functionality

- | | | |
|----|----------------------|--|
| 01 | AI symptom navigator | Collects symptoms and context, then suggests urgency and an appropriate care type with safety warnings. |
| 02 | Hospital directory | Searches verified facilities by location, hospital level, service, specialist, and operational availability. |
| 03 | Care matching | Connects assessment needs to suitable hospitals and departments instead of showing an unfiltered list. |
| 04 | Consultation request | Lets a first-time user verify contact details, submit care needs, and receive a trackable reference. |
| 05 | My Care | Keeps appointments, messages, medical history, notifications, and consent- |

controlled records in one place.

Scope and safety boundaries

The first version focuses on discovery, preliminary guidance, referrals, and consultation intake. It directs emergencies to 911 or the nearest emergency room, labels AI output as non-diagnostic, protects health data, and requires professional review for clinical decisions.

Project Proposal | 3

3. Mockup / Wireframe

User flow: Sign in > describe symptoms > review guidance > find matching care



CareNavigator
PH

Welcome back

Sign in to access your CareNavigator workspace.



[Forgot password?](#)

Sign in

[Create an account](#)

[Continue as guest](#)

Patient accounts are created or approved by authorized doctors. Hospital administrators create doctor accounts.

AI symptom navigator

Describe what you are experiencing. CareNavigator will guide you to the right level of care based on the urgency and the type of care to seek—it will not



Do not wait for this tool during an emergency. Go to the nearest emergency room now.

Tell us what you feel

Symptoms *

How long?

Existing conditions

Enter one per line or separate with commas.

Allergies

Enter one per line or separate with commas.

Current medications

Enter one per line or separate with commas.

Check symptoms



Home



Hospitals



Consult

SCREEN 1 - SIGN IN

Purpose: Provides secure entry for patients and care providers.

Key elements: Email, password, password recovery, and sign-in action.

Next action: Sign in and proceed to the symptom navigator.

Design notes

The interface uses a calm blue-and-teal palette, plain-language instructions, high-contrast actions, and persistent mobile navigation. Emergency messaging appears before the assessment to prevent unsafe reliance on the tool.

3. Mockup / Wireframe (continued)

for further evaluation and treatment

Recommended department

Primary Care

Warning signs

- Difficulty breathing
- Chest pain
- Severe headache

Possibilities to discuss with a clinician

- **Upper Respiratory Infection**
Fever, dry cough, and headache are common symptoms of upper respiratory infections
- **Influenza**
Fever, dry cough, and headache are common symptoms of influenza

Look for a hospital with

- None

Nearby hospitals for your needs

 **Level 1 or higher recommended: basic hospital and first-contact care.**

Location permission is blocked. Enable it in device settings to see nearby hospitals.

 [Try location again](#)

This is not a diagnosis, please consult a healthcare professional for proper evaluation and treatment

 [Find matching hospitals](#)

[Browse the full hospital directory](#)



Home



Hospitals



Consult



My care

Hospital director

Verified facilities, services, spe



Search hospital, loca



Show hospitals with availa

Hospital level:

✓ All

L



Level 3 • Tertiary Hospital

Tarlac Provincial H

 Tarlac City, Tarlac

 ER Available

 0 beds

[View hospital](#)



Level 2 • Secondary Hospital

President Ramon M Hospital

 Iba, Zambales

 ER limited

 0 beds



Home



Hospitals

Flow logic



Expected result

At the end of the flow, the user understands the suggested level of care and has a practical next step: view a suitable hospital, get directions, or begin a consultation. The design reduces uncertainty without removing professional judgment.

Wireframes 3-4 | 5

3. Mockup / Wireframe (continued)

Final step: begin a first-time online consultation



CareNavigator
PH



Home



Hospitals



Consult



My care



Emergency? Call 911 or
go to the nearest ER.

First-time online consultation

Verify your mobile number, tell us what you need, and receive a tra



Do not wait for an online consultation for severe breathing
to the nearest ER.

1

Verify your mobile number

Mobile number



+63

SMS delivery requires the S

SCREEN 5 - CONSULTATION REQUEST

Purpose: Gives the user a direct way to request professional care after finding an appropriate facility or department.

Key elements: Emergency warning, identity or contact verification, care-needs intake, schedule details, and a trackable consultation reference.

Next action: Verify the contact method, complete the intake form, and submit the consultation request.

Design rationale: This final screen closes the navigation loop. The app connects the user to a professional care workflow while continuing to display emergency safety instructions.

Complete five-screen user flow

1 SIGN IN > 2 ASSESS > 3 GUIDANCE > 4 HOSPITAL > 5 CONSULT

Wireframe 5 | 6

Appendix A. Complete Application Page Inventory

The five wireframes above satisfy the assignment's required user flow. The following appendix documents the remaining distinct mobile and web pages already implemented or represented in the CareNavigatorPH route structure.

Public access	Login, registration, home, symptom assessment, hospital directory, hospital detail, nearby map, and first-time consultation.
Patient access	Dashboard/My Care, care workspace, profile, secure messages, and notifications.
Doctor access	Clinical workspace, consultations, patients, result review, records, clinical actions, messages, and professional profile.
Hospital admin	Hospital profile, operations, consultations, hospital patients, doctors, schedules, availability, analytics, and alerts.
Super admin	Hospitals, administrators, service categories, doctors, reports, platform settings, permissions, audit, and security review.

Application routes reviewed

/login, /register, /home, /hospitals, /hospitals/map, /hospitals/:hospitalId, /assessment, /consult, /dashboard, /care, /profile, /messages/:conversationId, /notifications, /admin, and /admin/operations.

How to read the appendix

Each page below shows one distinct finalized interface or role-specific responsive view. Explanations identify the page's purpose, major functions, and intended user. Repeated test captures and duplicate screenshots are intentionally omitted.

A1. Web Login



CareNavigator
PH



Home



Hospitals



Consult



My care



Emergency? Call 911 or
go to the nearest ER.



CareNavigator
PH

Welcome

Sign in to access your account



Email



Password

Patient account
or Health Insurance ID

Audience: All registered users.

Purpose: Authenticates patients, doctors, hospital administrators, and super administrators before protected pages are opened.

Main functions: Email/password entry, password visibility control, password recovery, role-aware redirection, and access to account registration.

A2. Web Home



CareNavigator
PH



Home



Hospitals



Consult



My care



Emergency? Call 911 or
go to the nearest ER.

HEALTHCARE NAVIGATION FOR EVERY FILIPINO

Know where to go. Get care with confidence

Find verified hospitals, check live availability, and connect with a healthcare professional.

 Find a hospital

 Consult online



Severe breathing difficulty, chest pain, stroke signs, heavy

How can we help?



Find the right hospital

Compare services, specialists, ER status, beds, and rooms.

[Browse hospitals](#)



Your care in one place

Review consultations, records, results, and prescriptions securely.

[Open my care](#)



Consult online

Verify
online

[Start](#)

Audience: Guests and signed-in users.

Purpose: Serves as the application's navigation hub.

Main functions: Direct access to hospital search, online consultation, symptom assessment, My Care, and visible emergency instructions.

A3. Web Symptom Assessment



CareNavigator
PH



Home



Hospitals



Consult



My care



Emergency? Call 911 or
go to the nearest ER.

AI symptom navigator

Describe what you are experiencing. CareNavigator will suggest un



Do not wait for this tool during an emergency. Call 911 or

Tell us what you feel

Symptoms *

How long?

Age

Existing conditions

Allergies

Current medications

✦ Check symptoms

Audience: Guests or patients seeking preliminary guidance.

Purpose: Collects symptoms and health context in a structured form.

Main functions: Symptom description, duration, age, conditions, allergies, current medicines, emergency escalation, and assessment submission.

A4. Web Assessment Result

The screenshot shows the CareNavigator PH web assessment interface. On the left is a sidebar with a logo and navigation links: Home (highlighted), Hospitals, Consult, and My care. At the bottom of the sidebar is an emergency warning. The main area contains input fields for 'How long?' (2 days), 'Age' (30), 'Existing conditions', 'Allergies', and 'Current medications', each with a placeholder instruction. A large blue button at the bottom right is labeled 'Check symptoms'.

CareNavigator PH

- Home
- Hospitals
- Consult
- My care

Emergency? Call 911 or go to the nearest ER.

How long? 2 days

Age 30 30/4000

Existing conditions

Enter one per line or separate with commas.

Allergies

Enter one per line or separate with commas.

Current medications

Enter one per line or separate with commas.

Check symptoms

Audience: Users who completed the symptom form.

Purpose: Summarizes preliminary urgency and care-direction guidance.

Main functions: Suggested department, warning signs, possibilities to discuss with a clinician, non-diagnostic disclaimer, and hospital matching.

A5. Web Hospital Directory



CareNavigator
PH



Home



Hospitals



Consult



My care



Emergency? Call 911 or
go to the nearest ER.

Hospital directory

Verified facilities, services, specialists, and availability in one view.



Search hospital, location, service, or tag


Approved hospitals will ap

Audience: Anyone searching for care.

Purpose: Compares verified hospitals in one searchable directory.

Main functions: Search, level filters, ER availability, facility cards, locations, beds, rooms, services, and hospital-detail navigation.

A6. Web Hospital Detail



Audience: Users evaluating a specific hospital.

Purpose: Presents detailed information before travel or consultation.

Main functions: Hospital classification, ER/beds/rooms status, overview, location, contacts, departments, services, doctors, schedules, and directions.

A7. Web Nearby-Care Map



CareNavigator
PH



Home



Hospitals



Consult



My care



Emergency? Call 911 or
go to the nearest ER.

Nearby care



Use my location



ER available

Verified hospitals are ranked by straight-line distance. Travel time ca



Needed department, service, or specialty

Audience: Users choosing nearby care.

Purpose: Ranks suitable hospitals using the user's location and care requirement.

Main functions: Location permission, ER-only filter, required department/service entry, match ranking, and directions.

A8. Web Notifications



CareNavigator
PH



Home



Hospitals



Consult



My care



Emergency? Call 911 or
go to the nearest ER.

Audience: Signed-in patients, doctors, and administrators.

Purpose: Centralizes appointment, consultation, clinical, and operational alerts.

Main functions: Protected access, notification review, read status, and navigation to the related workflow.

A9. Doctor Patient Workspace



CareNavigator
PH



Overview



Clinical



Hospitals



Profile



My care workspace

Maria Santos · Doctor



Guest review



Consultations



Patients



Result review



Create patient



Juan Dela Cruz

CNPH-DEMO-0001 · Active

Hospital

Jose B. Lingad Memorial General Hospital

Blood type

O+

Conditions

Hypertension

Allergies

Penicillin



Care history across hospitals

Doctor-confirmed records shared with the assigned care team.



Respiratory consultation

Dec 19, 2025 · Bataan General Hospital and Medical Center

Attending physician: Dr. Elena Reyes

Diagnosis: Mild persistent asthma

Plan: Start controller inhaler and return for reassessment.



Hypertension follow-up

May 19, 2025 · Bataan General Hospital and Medical Center

Attending physician: Dr. Elena Reyes

Diagnosis: Essential hypertension

Plan: Continue maintenance medicine, home BP log, and low-sodium



Add record



Upload result



Prescription



Treatment plan



Medical document



Audience: Authorized doctors.

Purpose: Gives clinicians a consolidated view of assigned patients and cross-hospital care history.

Main functions: Patient creation, consultation history, diagnoses, treatment plans, prescriptions, laboratory requests, results, medical documents, and attachments.

A10. Doctor Profile



CareNavigator

PH



Overview



Clinical



Hospitals



Profile



My profile

Maria Santos · Doctor

Dr. M



Account information

Email	Doctor@demo.test
Mobile number	09170000003
Address	Not provided
Account status	Active



Professional profile

Specialization	Internal Medicine
Department	Internal Medicine
Hospital	Jose B. Lingad Memorial General Hospital
License number	DEMO-PHYSICIAN-001
Availability	Available
Consultation fee	₱500
Biography	Synthetic clinician profile for demonstrations

Audience: Doctors.

Purpose: Displays account and verified professional information.

Main functions: Contact information, account status, specialization, department, hospital assignment, license number, availability, consultation fee, and biography.

A11. Hospital-Admin Care Console

The screenshot displays the 'My care workspace' interface for a Hospital Admin. On the left is a sidebar for 'CareNavigator PH' with navigation links: Overview, Hospital, Operations, Care (highlighted), and Alerts. At the bottom of the sidebar is an emergency notice: 'Emergency? Call 911 or go to the nearest ER.' The main content area shows a consultation for 'Juan Dela Cruz' at 'Jose B. Lingad Memorial General Hospital' on 'Jul 21, 2026 6:00 PM'. The consultation description is 'Demo follow-up for blood pressure monitoring'. Below this, there are three buttons: 'Message', 'In Progress' (which is highlighted with a blue background), and 'Cancelled'.

**CareNavigator**
PH

 Overview

 Hospital

 Operations

 **Care**

 Alerts

 **Emergency? Call 911 or go to the nearest ER.**



**My care workspace**
Hospital Admin · Hospital Admin

 **Consultations**

 Hospital patients

Juan Dela Cruz

Jose B. Lingad Memorial General Hospital · Face To Face

Jul 21, 2026 6:00 PM

Demo follow-up for blood pressure monitoring

 **Message**

 **In Progress**

 **Cancelled**

Audience: Hospital administrators.

Purpose: Coordinates hospital-level consultation operations without unrestricted clinical-record access.

Main functions: Consultation monitoring, patient list, messages, status changes, hospital operations, availability, and alerts.

A12. Super-Admin Platform Console

 CareNavigator
PH

 Overview

 Platform

 Operations

 Alerts

 Emergency? Call 911 or go to the nearest ER.

  Administration
Demo Admin · Super Admin

 Hospitals

 Administrators

Hospital directory management

7 registered hospitals



Jose B. Lingad Memorial General Ho

Tertiary Hospital · City of San Fernando · Pampanga

Verified

Open



Bataan General Hospital and Medical

Tertiary Hospital · Balanga City · Bataan

Verified

Open



Bulacan Medical Center

Tertiary Hospital · Malolos City · Bulacan

Verified

Open



Dr. Paulino J. Garcia Memorial Resea

Tertiary Hospital · Cabanatuan City · Nueva Ecija

Verified


Open

Audience: Platform super administrators.


Purpose: Governs the CareNavigatorPH network.

Main functions: Hospital approval and management, administrator accounts, service categories, doctors, platform reports, permissions, settings, audits, and security review.

Appendix B. Role-Specific Mobile Views



B1. DOCTOR MOBILE CARE





B2. HOSPITAL-ADMIN MOBILE

Shows the doctor workspace in a compact layout, prioritizing clinical navigation, patient access, and mobile-safe actions.

Shows hospital operational access on a small screen, care coordination and administrative navigation.

Appendix B. Role-Specific Mobile Views (continued)



B3. SUPER-ADMIN MOBILE

Demonstrates the platform-management interface at mobile width and the responsive handling of administrative information.

B4. DESKTOP-PORTAL SAFEGUARD

Explains that selected high-risk administrative workflows require a desktop web portal, preventing an unsafe or unusable mobile experience.

4. Presentation Guide

Suggested 2-3 minute presentation script

OPENING	Our project is CareNavigatorPH, a mobile app that helps Filipino patients navigate from symptoms to suitable local care.
THE PROBLEM	People often do not know which hospital has the right service, specialist, or current capacity. This causes confusion, delays, and unnecessary travel.
SDG	The project supports SDG 3: Good Health and Well-Being, particularly Target 3.8 on access to quality essential healthcare services.
THE SOLUTION	Users can enter symptoms, receive preliminary and safety-focused guidance, see the recommended type of care, and find verified matching hospitals.
USER FLOW	The five wireframes show the complete journey: sign in, symptom assessment, care guidance, hospital selection, and consultation request.
SAFETY	CareNavigatorPH does not diagnose. Emergency cases are directed to 911 or the nearest ER, and clinical decisions remain with licensed professionals.
CLOSING	Our goal is simple: reduce the time and uncertainty between needing help and reaching the right care.

Success measures for a future prototype test

- Users can identify the next care step without assistance.
- Users can find a suitable facility in three minutes or less.
- Users understand that AI guidance is preliminary and non-diagnostic.
- Users notice emergency instructions before starting an assessment.

SUBMISSION CHECKLIST

- [] Replace the student, section, instructor, and date placeholders.
- [] Rename the file using your real LastName_FirstName.
- [] Review the wireframes and rehearse the presentation script.
- [] Upload the Word file or exported PDF as instructed.